

RNA QUANTIFICATION WITH CRYSTAL RT-DPCR

MULTIPLEXING TARGETS WITH CRYSTAL RT-DPCR

Reverse-transcription PCR (RT-PCR) is widely used for RNA analysis, for example for genomic expression or viral detection, and is applied to many fields, from life science research to diagnostics. Crystal RT-dPCR enables straightforward multiplexing of targets (references or transcripts of interest), absolute quan-

tification of viral RNA without need for standard curves and allows exquisite detection of fine differences in gene expression.

Thanks to 3-color Crystal RT-dPCR, three different messenger RNA (mRNA) can be analyzed in a single one-step reaction, demonstrating that 1.5-fold differences in gene expression can easily and precisely be quantified (**Figure 1**).




Figure 1: Absolute quantification of mRNAs in total RNA using one-step Crystal RT-dPCR. Glucuronidase beta (GUSB), Albumin (ALB) and Translin (TSN) mRNAs were quantified by 3-color Crystal RT-dPCR in a serial dilution of total human RNA, ranging from 0.06 to 3.6 ng/ μ L with a ratio of 1.5 between each point. Assays were performed in Sapphire Chips using 1x qScript™ XLT One-Step RT-qPCR Toughmix®, and FAM-, HEX- and Cy5-labeled TaqMan™ probes. Results for the three replicates show good linearity, accuracy and reproducibility for the different levels of expression of each target. Slope of the regression line: GUSB = 588.8 cp/ng; ALB = 96.5 cp/ng; TSN = 980.6 cp/ng.

QUANTIFICATION OF LOW-LEVEL TRANSCRIPTS

Detecting low amounts of transcripts in high background of total RNA can be challenging. Quantification with the naica® sys-

tem shows excellent linearity and reproducibility for low fraction of transcript (**Figure 2A**). Reaction efficiency, as assessed through fluorescence amplitude, remains stable even in presence of a high load of total RNA (2 μ g - **Figure 2B**).




Figure 2: Absolute quantification of mRNA with a high amount of total RNA using one-step Crystal RT-dPCR. Albumin (ALB) mRNA was quantified by Crystal RT-dPCR in a serial dilution of total human RNA, ranging from 2 μ g to 125 ng per reaction with a ratio of 2 between each point. Assays were performed in Sapphire Chips using 1x qScript™ XLT One-Step RT-qPCR Toughmix® and HEX-labeled TaqMan™ probe. Linear plot shows accurate and reproducible results for the three duplicates (**A**), and no loss of efficiency has been observed on droplets fluorescence dotplots (**B**). RFU: Relative Fluorescence Units. *27 μ L per reaction.

VIRAL RNA ABSOLUTE QUANTIFICATION

Dengue virus 1 (DEN-1) and its 3 other closely-related serotypes can cause fatal conditions, such as Dengue Hemorrhagic Fever and Dengue Shock Syndrome. No need for calibration curves using standards: in less than 2h30min, Crystal Digital RT-PCR enables viral quantification in a fast and user-friendly manner, suitable to all labs and workflows.

To learn more about digital PCR, please visit Stilla Technologies' Learning Center at stillatechnologies.com/digital-pcr




Figure 3: Absolute quantification of Dengue-1 RNA using one-step Crystal RT-dPCR. Dengue-1 RNA was quantified by Crystal RT-dPCR in a background of 300 cp of total human RNA, with ALB mRNA quantified as a reference. Assays were performed in Sapphire Chips using 1x qScript™ XLT One-Step RT-qPCR Toughmix® and FAM- and HEX-labeled TaqMan™ probes. Dengue-1 RNA was previously quantified by the manufacturer by digital PCR. Logarithmic plots show good concordance with the expected concentration and reproducible results.